



BIOTECH
OPHTHALMICS

VITREORETINAL

A NEW GENERATION OF
VITREO-RETINAL TAMPONADES

Bio Decalin

DESCRIPTION

- + Bio Decalin is sterile, non-pyrogenic ultra-purified 100% fluorinated Perfluorodecalin liquid as an aid for retinal surgery. It is chemically inert and non-toxic in nature.

INDICATION

Intra-operative device for retinal surgery. It is used in the surgical management of:

- + Retinal Detachment with Giant Tears
- + Retinal Detachment with Proliferative Vitreoretinopathy (PVR) including Diabetic Retinopathy
- + Retrieval of luxated natural or artificial lens in the Vitreous.

ADVANTAGES

- + Fewer bubbles during surgery because boiling point is 142°C
- + Effective tamponade with Silicone Oil in complicated cases

SUPPLY

Bio Decalin is supplied in a glass vial as well as pre-filled syringe & Vial delivering 5.0 ml / 7.0 ml of sterile Perfluorodecalin liquid. The product is packed in medical grade blister packaging and finally supplied along with information for use (IFU).



BIOTECH
OPHTHALMICS

Factory

1st Floor, 555 / 556 / 557 , Nr. Subham Tex-O-Pack,
Khatraj, Dist. Gandhinagar, Gujarat, India.

India Business

210, Sarthik - II, Opp, Rajpath Club,
Sarkhej-Gandhinagar Highway,
Ahmedabad-380054, Gujarat, India.
Ph. / Fax: 079-40038374,32420133
E-mail: sales@biotechvisioncare.com

International Business

401, Sarthik- II, Opp. Rajpath Club,
Sarkhej-Gandhinagar Highway,
Ahmedabad-380054, Gujarat, India.
Ph.: +91-79-26870896/97/98 **Fax:** +91-79-40033767
E-mail: intbusiness@biotechvisioncare.com

Web: www.biotechvisioncare.com