

CE  
2460

# VITREORETINAL

A NEW GENERATION OF VITREO-RETINAL TAMPONADES

biotech

VISIONCARE

## Bio Blue 90

An adjunct to Vitreo-retinal Surgery

---

### DESCRIPTION

- + Bio Blue 90 is a sterile, pyrogen-free, non toxic, isotonic and ready to use solution of Acid Blue 90 or Brilliant Blue G (a non-Fluorescent Triphenylmethane stain) with effective staining properties.

### INDICATION

- + Ready to use sterile solution - no mixing or filtration is required.
- + Excellent biocompatibility - a more favorable safety profile.
- + Better affinity to ILM and low affinity to ERM and hence better visualization of ILM leading to ease of peeling during macular hole and macular pucker surgery.
- + Facilitate the successful creation of CCC (Continuous Curvilinear Capsulorhexis).
- + No adverse / sideeffect. Better retinal tolerance.
- + No side reaction such as photo induced cross linking of collagen fibres.
- + Easy to remove the residual product from the eye.

### SUPPLY

Supplied in a 0.5 ml pre-filled syringe for easy opening and use.



**BIOTECH VISION CARE PVT LTD** | Block 1, Abhishree Corporate Park, Opp. Swagat Bungalow BRTS Stop  
Bopal - Ambli Road, Ahmedabad - 380 058 Gujarat India. | Tel: +91 79 66823000  
Email: [indiasales@biotechhealthcare.com](mailto:indiasales@biotechhealthcare.com) | [intlsales@biotechhealthcare.com](mailto:intlsales@biotechhealthcare.com) | [www.biotechhealthcare.com](http://www.biotechhealthcare.com)