



BIOTECH
OPHTHALMICS

VITREORETINAL

A NEW GENERATION OF
VITREO-RETINAL TAMPONADES

Biosil / Biosil -F

A reliable and proven treatment for complex Retinal Detachments

DESCRIPTION

- + Biosil / Biosil-F is an ultra-purified Silicone Oil which delivers a maximum interfacial tension and minimizes interactions between tissues, cells and endo-tamponades media. The physical properties include a combination of specific gravity, refractive index and surface tension.

ADVANTAGES

- + The choice of viscosity offers an optimum balance between easy injection and a stable temporary tamponade
- + Less tissue impregnation
- + Highly purified and low level of reactive - OH end groups results in reduced incidence of emulsification
- + Non-pyrogenic, clear and colorless
- + Relatively inert and preservative free
- + Minimal volatile contents

PRIMARY INDICATIONS

- + Proliferative Vitreoretinopathy (PVR)
- + Retinal Detachment caused by trauma
- + Giant Retinal Tears (GRTs)
- + Diabetic Tractional Detachment
- + Cytomegalovirus (CMV) Retinitis

Properties	Biosil	Biosil-F
Viscosity	1000 cSt	5000 cSt
Specific Gravity	0.972 g/cm at 25° C	0.972 g/cm at 25° C
Refractive Index	1.40	1.40
Packing	10.0 ml Vial / 10.0 ml PFS	10.0 ml Vial / 10.0 ml PFS



BIOTECH
OPHTHALMICS

Factory

1st Floor, 555 / 556 / 557 , Nr. Subham Tex-O-Pack,
Khatraj, Dist. Gandhinagar, Gujarat, India.

India Business

210, Sarthik - II, Opp, Rajpath Club,
Sarkhej-Gandhinagar Highway,
Ahmedabad-380054, Gujarat, India.
Ph. / Fax: 079-40038374,32420133
E-mail: sales@biotechvisioncare.com

International Business

401, Sarthik- II, Opp. Rajpath Club,
Sarkhej-Gandhinagar Highway,
Ahmedabad-380054, Gujarat, India.
Ph.: +91-79-26870896/97/98 **Fax:** +91-79-40033767
E-mail: intbusiness@biotechvisioncare.com

Web: www.biotechvisioncare.com